




NS400BRUSH is special wand, with a cylindrical rotating brush, to be used for cleaning carpet in combination with an injection-extraction machine. The mechanical brushing action allows to remove the dirt from the base of the carpet fibres. No residual dirt will re-emerge on the carpet surface a few days after cleaning. Even a very dirty carpet is left completely clean in a single pass. It is not necessary to use a lot of water to clean carpet, therefore the drying time is minimal: carpet can be dry in 1 hour.



Injection-extraction cleaning with rotating brush



Advantages of NS400BRUSH

1. The result of injection-extraction cleaning with rotating brush is excellent.
2. Since the dirt is completely removed from the base of carpet fibres, no residual dirt will re-emerge on carpet surface after few days, as it happens with surface-cleaning systems!
3. Carpet is left almost dry after cleaning and can be completely dry in 1 hour.
4. A lever on the handle allows to adjust the water-flow, for a deep or a maintenance cleaning.
5. The "S-shape" allows to easily clean as well under tables and desks.
6. It is possible to use the wand NS400BRUSH with any Santoemma injection-extraction machine.

Technical Data

Cleaning width	400 mm
Number of nozzles:	4
Brush:	120 W 600 rot/min 364 mm
Cable length:	15 m
Weight	16 Kg
Weight (packaged wand):	18 Kg
Dimensions:	cm 42 x 42 x 90 (h)
Package dimensions:	cm 48 x 65 x 47 (h)

Complete wand NS400BRUSH, for injection-extraction cleaning with rotating brush.